

## Disponibilité des analytes, programme QMEQAS

Analyte		QMEQAS			
		Sang	Urine	Sérum	Cheveux
Aluminium	Al	Inclus	Inclus	Inclus	Inclus
Antimoine	Sb	Inclus	Inclus	Inclus	Inclus
Argent	Ag	Inclus	Inclus	Inclus	Inclus
Arsenic total	As	Inclus	Inclus <sup>1</sup>	Inclus	Inclus
Baryum	Ba	Inclus	Inclus	Inclus	Inclus
Béryllium	Be	Inclus	Inclus	Inclus	Inclus
Bismuth	Bi	Inclus	Inclus	Inclus	Inclus
Bore	B		Inclus		
Cadmium	Cd	Inclus	Inclus	Inclus	Inclus
Césium	Cs	Inclus	Inclus		
Chrome	Cr	Inclus	Inclus	Inclus	Inclus
Cobalt	Co	Inclus	Inclus	Inclus	Inclus
Cuivre	Cu	Inclus	Inclus	Inclus	Inclus
Étain	Sn	Inclus	Inclus	Inclus	Inclus
Iodure	I	Inclus	Inclus	Inclus	Inclus
Lithium	Li		Inclus		
Manganèse	Mn	Inclus	Inclus	Inclus	Inclus
Mercure	Hg	Inclus	Inclus	Inclus	Inclus
Molybdène	Mo	Inclus	Inclus	Inclus	Inclus
Nickel	Ni	Inclus	Inclus	Inclus	Inclus
Platine	Pt	Inclus	Inclus	Inclus	Inclus
Plomb	Pb	Inclus	Inclus	Inclus	Inclus
Sélénium	Se	Inclus	Inclus	Inclus	Inclus
Strontium	Sr	Inclus	Inclus		
Tellure	Te	Inclus	Inclus	Inclus	Inclus
Thallium	Tl	Inclus	Inclus	Inclus	Inclus
Thorium	Th	Inclus	Inclus	Inclus	Inclus
Titane	Ti	Inclus	Inclus	Inclus	
Uranium	U	Inclus	Inclus	Inclus	Inclus
Vanadium	V	Inclus	Inclus	Inclus	Inclus
Zinc	Zn	Inclus	Inclus	Inclus	Inclus

1. La spéciation de l'arsenic est disponible dans l'urine. Sept formes sont disponibles, soit total, inorganique, DMAA, MAA, As+3, As+5 et alimentaire.